

# BioPAK Computer Requirements

These are the minimum requirements to operate the BioPAK Web v1.0.0 software. Please ensure your computer meets these minimum requirements before purchasing the BioJVA Wireless device.

Specification	With Adapter	Without Adapter
Operating System	Windows 10 or later	
	MacOS OS X 10.10 Yosemite or later*	
	ChromeOS (any version)	
	Android 6.0 Marshmallow or later	
Computer Browser	Google Chrome version 78 or later	
	Microsoft Edge version 79 or later	
	Opera version 65 or later	
	Samsung Internet version 12 or later	
	Android Chrome version 78 or later	
	Android Opera version 56 or later	
Link Manager Protocol (LMP)	4.0	9.0
Memory (RAM)	4 GB	4 GB
Disk Storage	20 GB available	20 GB available
Processor	Intel Pentium 4, AMD Athlon 64, or newer	

\*The wireless adapter is only compatible with Windows 10 or later.